

© 2015 HORIBA, Ltd. All rights reserved.

HORIBA Scientific Patrick CHAPON

Direct measurement of GDOES crater depth with built-in Differential Interferometer

2015 HORIBA, Ltd. All rights reserve

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific





May 2017

HORIBA



HORIBA Scientific: European Operation

Come to visit !



HORIBA Paris Saclay Site

© 2015 HORIBA, Ltd. All rights reserved.

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific



Explore the future



Thickness? Spectroscopic Ellipsometry





Output of a GD-OES analysis



R. Escobar Galindo, E. Forniès and J. M. Albella (2005). *J. Anal. At. Spectrom.*, Vol 20, pp 1131-1138

2015 HORIBA, LIG. All rights reserved.		
Explore the future	Automotive Test Systems Process & Environmental Medical Semiconductor Scientific	HORIBA



Principle of GD-OES



- In analytical GDs the sample is mounted outside the plasma chamber.
 - The Ar plasma is confined inside a cylindrical anode and the sample is placed outside, acting as the cathode.
- Ar+ ions are accelerated towards the surface of the sample and they sputter it.
- The same plasma excites the sputtered species and light emission is observed.
- The real time measurement of the emitted light with an optical spectrometer therefore gives a depth resolved profile of the investigated material.

© 2015 HORIBA, Ltd. All rights reserved. Explore the future 6



Until now depth = calculation



Pure Ti layers deposited by two different sputtering techniques \rightarrow different density – not easily taken in account by software.

© 2015 HORIBA, Ltd. All rights reserved.

Explore the future





The Analytical Scientist – May 2014



What are the opportunities and challenges?

GH: There are plenty of challenges (or opportunities, depending on your point of view):

- Low-noise, wide-bandwidth single-photon detectors
- · Tunable lasers for the entire UV-Vis wavelength range
- Broad-scene imaging at nanometer spatial resolution
- Bright non-laser sources for IR
- · Low noise, ultra-wide bandwidth analog amplifiers
- The need to develop a method for measuring sample erosion on the nanometer scale.

Gary Hieftje, Peter Griffiths, and Volker Deckert share their views on the past, present, and bright future of spectroscopy

© 2015 HORIBA, Ltd. All rights reserved.

Explore the future

8



Constraints (2) and Challenges (2)



Technical solution: DIP

- Interferometric method.
- Relative measurement between the crater and the surface close to the crater.
- Red laser diode (~635 nm +/- 5 nm), no interesting spectroscopic line in this spectral region.







Outputs from the interferometer



2015 HORIBA, Ltd. All rights reserved.		11/16
Explore the future	Automotive Test Systems Process & Environmental Medical Semiconductor Scientific	HORIBA



Experimental setup



The detection scheme shown in Figure 3 is that of a homodyne polarization interferometer [12]. After the reflection onto the sample, the two laser beams are recombined. The different polarizations (linear and circular) are then separated and measured by using four photodiodes (D1 to D4, in Figure 3), giving signals 51 to 54 of the form

$S_1 = A_1 + B_1 \cos \Phi$
$S_2 = A_2 - B_2 \cos \Phi$
$S_3 = A_3 + B_3 \sin \Phi$
$S_4 = A_4 - B_4 sin \Phi$

where the parameters A_i and B_i may be slowly varying functions of time.

For non transparent samples, the crater depth D is linked to the phase ϕ according to

 $D = \frac{1}{4\pi} \Phi$ Eq.(1)

Ref: article submitted in JAAS

Figure 3. Schematic representation of the optical path for the DiP setup, viewed from the top (a), and the side (b). L: laser source, W: Wollaston polarization beam splitter, R: 45° polarization rotator, (N)PBS : (non) polarizing beamsplitters, QWP : quarter waveplate, D1 to D4 : Silicon photodiodes.

© 2015 HORIBA, Ltd. All rights reserved.		
Explore the future	Automotive Test Systems Process & Environmental Medical Semiconductor Scientific	HORIBA



Case 1: Bulk sample

The simplest case: Only one reflecting surface





Figure 7. GD-OES analysis of a Si wafer. The Ar gas pressure was kept constant and the power was randomly varied from 20 W to 70 W. (a) Evolution of the Si signal as a function of the sputtering time. It shows the rapid response and the variation of the light intensity when the RF power is varied. (b) The evolution of the desth as a function of the sputtering time. For each variation of the applied power, a variation in the desth solve is observab. Errosion rates are reported in the figure.

© 2015 HORIBA, Ltd. All rights reserved.

Explore the future





Case 2: Multilayered sample

Only one reflecting surface at a time



$$\Delta \varphi(t) = \Delta \varphi^{depth}$$

$$\int$$

$$\Delta \varphi^{depth} = \frac{4\pi}{\lambda} \times d$$
Only depth related
$$\Rightarrow d = \frac{\Delta \varphi}{4\pi} \times \lambda$$

$$\Rightarrow \text{ conversion factor : 50.5 nm/rad}$$

© 2015 HORIBA, Ltd. All rights reserved. Explore the future





Example : PVD coating



© 2015 HORIBA, Ltd. All rights reserved.

Explore the future





Accuracy



Figure 6. Depth measured by DiP as a function of the depth measured by the mechanical profilometer. Solid line is a linear fit of equation y=1.003x-2.2.





Comparison with Ellipsometry

• GD + DIP

Ellipsometry

- Sputtering
- Fast Elemental Depth Profile
- Depth by Interferometry
- One laser wavelength
- Needs reflective surface/interfaces
- Best nm !

- Change of polarization state
- Very sensitive and ideal for complex materials (anisotropy, gradients etc); importance of modelling.
- Multi wavelengths
- Needs transparent/semi transparent layers
- Non destructive, possible in situ
- Best Å !

© 2015 HORIBA, Ltd. All rights reserved.



DIP Summary

Optical technique coupled to standard GD Differential Interferometry

Ideal for reflective samples

A new tool for thickness determination

© 2015 HORIBA, Ltd. All rights reserved.

Explore the future

Automotive Test Systems | Process & Environmental | Medical | Semiconductor | Scientific

HORIBA

18



Thank you very much for your attention

© 2015 HORIBA, Ltd. All rights reserved.

Explore the future

