

ThermoFisher SCIENTIFIC

Calibration factors in Fast Flow Glow Discharge Mass Spectrometry (FF-GD-MS): continuous versus pulsed mode

Joachim Hinrichs 17. GDS-Anwendertreffen, Bremen, May 2017 Thermo Fisher Scientific, Bremen, Germany

The world leader in serving science

Scope

- Many GD-MS analyses rely on semiquantitative analyes
- Based on a general table of Relative Sensitivity Factors
- Universally applied to all matrices = Standard RSF
- Results of some elements tend to be low – mostly Ti and V
- Can pulse mode use the same RSF factors than continuous DC mode?



 Note: for many materials, CRM can be used. This yields accuracies close to solution based calibration by e.g. ICP-MS.



Poster at the Plasma Winter Conference 2017

thermo scientific

Calibration factors in Fast Flow Glow Discharge Mass Spectrometry (FF-GD-MS): continuous versus pulsed mode

Joachim Hinrichs, Shona McSheehy Ducos, Thermo Fisher Scientific, Bremen, Germany. joachim.hinrichs@thermofisher.com

ABSTRACT

Purpose: To compare and revise calibration factors for the Thermo Scientific[™] Element GD[™] Plus GD-MS in continuous mode and pulsed mode operation.

Methods: The Element GD Plus GD-MS was used to measure a set of 23 iron and Nickel based certified reference materials, selected to cover a wide range of elemental concentrations.

Results: Calibration factors for 33 elements were obtained from the set of reference materials analyzed. In continuous mode operation, the calibration factors confirm those previously used for semiguanitative analyses. Significantly improved activity (SANS expected for T), V, Cr, and Zr with the data obtained during this study. For plused mode, the linearity at high concentrations is improved. Nost elements are analyzed within the accuracy range of ±30% considered typical for semiguantitative GD-MS analyses. For improved accuracy, the USE Or a declarate set of calibration factors is advisable.

INTRODUCTION

Sector field glow discharge mass spectrometry is applied for the analysis of high purty bulk metals and alloys, semiconductors and ceramics, especially by the aerospace, electronics and photovoltaic industries. The Element GD Plus GD-MS features a fast flow glow discharge source that can be operated in continuous or pulsed mode. Continuous mode operation dires the advantage of high sputter rates to remove contaminated



Figure 1. Analytical procedure.



Table 2. Instrumental parameters (standard conditions).

Parameter	Continuous DC mode	Pulsed DC mode
GD Source Pressure [mbar]	2.8	2.8
Discharge Voltage (V)	-700	1000 (set value)
Discharge Current [mA]	25 (set value)	-11
Discharge Gas (mL/min)	-240	-240
Extraction [V]	-2000	-2000
Focus [V]	-1200	-1200
Pulse Duration (pSec)	0.8	50
Pulse Mode	ne	yes
Pulse Frequency [kHz]	\$\$	2

Figure 3 shows the calibrations for V and Cr, svalable at higher levels, with C neaching + 20% in a few mailwrias. At such high concentrations, matrix effects can occur, and a distinct difference between the measurement modes can be observed.

In continuous mode, V and especially Cr show a significant offset from the regression line in the low percent concentration range, while pulsed mode shows a much better fit. Since pulsed mode gives the same 835 over the entire range, this is therefore the preferable option for general survey analyses.

Figure 3. Calibration ourves for V and Cr at the percent level. Left: continuous mode, right: pulsed mode. Note the significantly better fit for pulsed mode analysis.



For pulsed mode analysis, only a small data base was so far available ⁴). Table 2 shows that the pulsed mode RSF overlap for many elements with continuous mode. In datail though, a number of elements show significantly lower RSF values, especially Mg, Ca, Mn, Zn, Se, Ag, Sn, 'Pe, Pband BL Others, e.g. SI, NI, As and So, yielb higher RSF. With these differences, the use of an RSF table derived from continuous mode operation may result in larger errors. Therefore, a dedicated RSF table for pulsed mode is recommended to support accurate routine operation for sembuantative analyses.

The standard instrumental parameters used for pulsed mode cover the vest majority of applications. Still, further studies to investigate the influence of other pulse parameters, e.g. pulse duration and frequency, need to be carried out. Initial results (not shown) indicate a stable instrumental response br pulse durations 2 40 µs and pulse frequencies in the range 1 to 4 kHz.

Figure 4. Calibration ourves for selected non-metallic elements. Left: continuous mode, right: pulsed mode.







ThermoFisher SCIENTIFIC Quantification by GD-MS

The world leader in serving science

Quantification by GD-MS

- Main application is concentration determination
- Concentration = mass fraction
- · Isotopes are measured, concentrations are calculated and reported
- Abundance correction for isotopes is automatically applied
- No blanks, GD considered as blank free
- IBR = Ion Beam Ratio. This means normalizing all signals to matrix signals.
 Acts like an Internal Standard in ICP. IBR used for RSF quantification
- RSF concept = Relative Sensitivity Factor. Evaluation for semiquant = recalculating calibration factors based on measured matrix composition
- A general RSF table is commonly used for semiquantitative analysis by GD-MS. This is not absolutely accurate, but at least consistent. This is of major importance for quality control in industrial labs.
- Calibrated RSF give a matrix specific calibration = best accuracy



Evaluation by GD-MS

Evaluation Glow Discharge

Basic Assumption

$$\sum_{Elements} c_y = 1$$

 $IBR_{x} = \frac{I_{x}}{\sum I_{y}}$

Elements

I = intensity measured [cps]

c = concentration

- *IBR* = ion beam ratio
- *RSF* = relative sensitivity factor.

 $(IBR_x \text{ all intensities known})$

$$RSF_x = \frac{c_x}{IBR_x}$$

 $\sum IBR_y = 1$

Elements

$$c_x = RSF_x \cdot IBR_x$$

$$c_x = \frac{RSF_x \cdot IBR_x}{\sum_{Elements} RSF_y \cdot IBR_y}$$



The RSF concept works because in GD-MS...

- ...the ionization is mostly occurring when the cloud of atoms is travelling in the Argon gas to the mass
- The GD plasma is mostly Argon with sample atoms, causing a very similar ionization regime even for different materials = similar environment for different materials
- The key is to consider the ionization efficiencies relative to each other, and maintain the same ionization process, which means similar Glow Discharge conditions
- For FF-GD-MS, the **pulsed** discharge generates a better equilibrated aerosol with uniform ionization. This is reflected by low influence of source parameters onto the IBR and Concentrations measured



Quantification by GD-MS

- Standard RSF are not fixed, but are instrument dependent: Source, Mass response of the MS, Tuning.
- In continuous DC mode, RSF were dependent on gas flow needs to be strongly controlled
- In **pulsed** mode, there is just a small dependency on gas flow
- Semiquant = STD RSF typically in a range of $\pm 30\%$ accuracy.
- E.g. 10 ppm as a semiquant result should be regarded as "between 7 and 13 ppm"
- For high purity materials, usually semiquant results are sufficient, e.g. for "<0.01 ppm"



The existing Standard RSF file lists sensitivity factors measured from steel and Iron matrix, normalized to 100% Fe:

Help _ & ×

2

					_			•					
📑 Results	- Thermo E	LEMENT -	[Standard	RSF.rsf			x		Results	- Thermo E	LEMENT	- [Standard	I RSF.rsf
🛐 <u>F</u> ile	<u>R</u> SF Table	<u>Options</u>	<u>V</u> iew <u>W</u>	(indow	<u>H</u> elp	_ d	s ×		🔝 <u>F</u> ile	<u>R</u> SF Table	<u>Options</u>	<u>V</u> iew <u>V</u>	<u>N</u> indow
🗃 🖬			المنا الم		?				<i> </i> 🖬			1.8 1.8	
RSF							*		Rb	1.0000	1.0000	1.0000	
									Sr	0.5200	0.5200	0.5200	
Element	LR	MR	HR						Y	0.5400	0.5400	0.5400	
									Zr	0.5600	0.5600	0.5600	
н	1.0000	1.0000	1.0000						Nb	0.6600	0.6600	0.6600	
Не	1.0000	1.0000	1.0000						Mo	0.9200	0.9200	0.9200	
Li	3.1000	3.1000	3.1000						IC	1.0000	1.0000	1.0000	
Be	5.0400	5.0400	5.0400						RU	0.7600	0.7600	0.7600	
В	6.4900	6.4900	6.4900				Ξ		Rn Dd	1.3600	1.3600	1.3600	
С	9.2700	9.2700	9.2700						Pa	1.8100	1.8100	1.8100	
N	1.0000	1.0000	1.0000						Ag	3.8500	3.8500	3.8500	
0	1.0000	1.0000	1.0000							3.4100	3.4100	3.4100	
F	1.0000	1.0000	1.0000						50	1.4500	1.4500	1.4000	
Ne	1.0000	1.0000	1.0000						500 Sh	1.2900	1.2900	1.2900	
Na	0.9500	0.9500	0.9500						30 To	4.0900	4.0900	4.0900	
Ma	1 5100	1 5100	1 5100							4.4300	4.4300	4.4300	
	1 2700	1 2700	1 2700						v.	1.0000	1.0000	1.0000	
AI Ci	3.0400	3.0400	3.0400						Ce Ce	1.0000	1.0000	1.0000	
31	2 6600	2 6600	2 6600							0.5600	0.5600	0.5600	
P c	2 4200	2 4200	2 4200						La	0.5000	0.5000	0.5000	
3	4,0000	3.4300	3.4300						Ce	0.7500	0.7500	0.7500	
	1.0000	1.0000	1.0000						Pr	0.8100	0.8100	0.8100	
	1.0000	1.0000	1.0000						Nd	0.8800	0.8800	0.8800	
<u>к</u>	0.5900	0.5900	0.5900						Pm	1.0000	1.0000	1.0000	
Ca	0.4500	0.4500	0.4500						Sm	0.7400	0.7400	0.7400	
SC	0.5900	0.5900	0.5900						Eu	0.6500	0.6500	0.6500	
	0.4100	0.4100	0.4100						Gd	0.6800	0.6800	0.6800	
V	0.5400	0.5400	0.5400						Tb	0.8300	0.8300	0.8300	
Cr	1.2800	1.2800	1.2800						Dy	0.7800	0.7800	0.7800	
Mn	1.0100	1.0100	1.0100						Ho	0.8100	0.8100	0.8100	
Fe	1.0000	1.0000	1.0000						Er	0.7600	0.7600	0.7600	
Co	1.0400	1.0400	1.0400						Tm	0.7500	0.7500	0.7500	
Ni	1.5100	1.5100	1.5100						Yb	0.8000	0.8000	0.8000	
Cu	2.4400	2.4400	2.4400						Lu	0.6700	0.6700	0.6700	
Zn	3.8300	3.8300	3.8300						Hf	0.7900	0.7900	0.7900	
Ga	2.3400	2.3400	2.3400						Та	1.2400	1.2400	1.2400	
Ge	1.3300	1.3300	1.3300						w	1.6100	1.6100	1.6100	
As	5.1300	5.1300	5.1300						Re	0.9700	0.9700	0.9700	
Se	3.7700	3.7700	3.7700						Os	1.8200	1.8200	1.8200	
Br	1.0000	1.0000	1.0000						Ir	1.3800	1.3800	1.3800	
Kr	1.0000	1.0000	1.0000				-		Pt	2.4300	2.4300	2.4300	
DK.	4 0000	4 0000	4 0000				-		Au	2.3700	2.3700	2.3700	
Ready									Ready				

- Continuous DC mode
- Discharge current 40
 45 mA
- NIST and BAM standards
- Full procedure was
 not published



Why revisiting the Standard RSF?

- Reports on problematic elements. Liverpool EW-GDS meeting 2016
- Gas flow dependency difficult to define best conditions
- What is the effect of the pulsed mode on RSF?

Revisions:

- 2004 original version based on BAM doped pellets plus some steel CRM NIST 176x series
- 2008 update for Carbon: RSF for C in steel ~ 9.3



Standard RSF file

Anal Bioanal Chem DOI 10.1007/s00216-006-0645-5

SPECIAL ISSUE PAPER

Ralf Matschat · Joachim Hinrichs · Heinrich Kipphardt

Application of glow discharge mass spectrometry to multielement ultra-trace determination in ultrahigh-purity copper and iron: a calibration approach achieving quantification and traceability

	BAM_Dot 3-6	BAM_Dot 7-9	Fe CRM	Average				
Ti		0.44	0.39	0.415				
V		0.59	0.49	0.540				
Cr	1.49	1.09	1.26	1.280				
Fe	1.00	1.00	1.00	1.00				
Zr		0.57	0.56	0.565				
Sn		1.19	1.4	1.295				
Sb		4.88	4.9	4.890				

Example.







Standard RSF re-visited

	Туре
ECRM098-1	High purity iron
ECRM270-1	High alloy steel 1.4835
ECRM271-1	Tool steel 1.2344
ECRM289-1	High temperature steel
ECRM295-1	Highly alloyed steel
ECRM298-1	Duplex stainless steel
ECRM297-1	Radionox steel 1.4696
ECRM379-1	Highly alloyed steel
NIST1261	AISI 4340 steel
NIST1262	AISI 94B17 steel
NIST1263	Cr-V steel
NIST1264	High carbon steel
NIST1265	Electrolytic iron
NIST1761	Low alloy steel
NIST1762	Low alloy steel
NIST1766	Low alloy steel
NIST1767	Low alloy steel
NIST1173	Ni-Cr-Mo-V steel
MBH12X353	Low alloy steel (wrought)
BS2205	Duplex alloy 2205
NIST1249	Nickel superalloy 718
IARM59C	Nickel alloy 825
BAS346A	Nickel alloy IN 100

Procedure:





Pulsed mode



EN

Example calibrations @ trace level

- 0 **X** rmo ELEMENT - [Fe JH161101 DC H.cal] Results - Thermo ELEMENT - [Fe_JH161101 p_H.cal Eile Calibration Display Options View Window Help @ Eile Calibration Display Options View Window Hel 📽 🖬 🖉 < 🗄 🖌 🐚 📆 🖉 😭 😤 🖬 🖉 < 🕒 🖌 🐚 🖏 🖉 🖓 😭 Ti47(MR) y = 0.56 * x Regr. Type : Thru Zero RSF : 0.5599 Ti47(MR) y = 0.38 * x Regr. Type : Thru Zero RSF : 0.3836 Continuous DC 47**T**i 47**T**i 1500 21 Ion Beam Ratio [ppm] 2000 3000 Ion Beam Ratio [ppm] Zr91(MR) y = 0.72 * x Regr. Type : Thru Zero RSF : 0.7180 Zr91(MR) y = 0.47 * x Regr. Type : Thru Zero RSF : 0.466 ⁹¹Zr ⁹¹Zr 1500 3000 Ion Beam Ratio (ppm Ratio Ir NUM Sb121(MR) y = 4.79 * x Regr. Type : Thru Zero RSF : 4.7924 Sb121(MR) y = 6.42 * x Regr. Type : Thru Zero RSF : 6.4245 ¹²¹Sb ¹²¹Sb 30 Ion Beam Ratio (ppm 40 m Ratio (pp NUM

Right graphs: **Pulsed Mode**



Left graphs:

Continuous DC Pulsed Mode



Less variation, better fit on a linear curve for especially Cr and V at higher level in pulsed mode.



Results: Selection of RSF

	Standard RSF	Continuous Mode RSF	Pulsed Mode RSF	# Cert. values	# Info values
В	6.49	8.66	6.40	16	-
С	9.27	11.2	7.22	22	-
Mg	1.51	1.49	0.80	11	5
Si	3.04	4.06	4.96	21	-
Р	3.66	4.69	4.15	20	-
Ti	0.41	0.56	0.38	16	3
V	0.54	0.76	0.51	21	-
Cr	1.28	1.53	1.20	22	-
Mn	1.01	0.96	0.57	22	-
Zr	0.56	0.72	0.47	11	4
Мо	0.92	0.89	0.75	22	-
Ag	3.85	3.25	2.00	6	3
Sn	1.29	1.24	0.71	18	3
Sb	4.89	4.79	6.42	16	5
Pb	1.36	1.28	0.86	13	5
Bi	2.94	2.59	2.00	7	2



Conclusion

- Step 1 Ti, V, Cr, Zr give higher Standard RSF values during these experiments.
- Step 2 verify at several instruments.
- Step 3 generate new revision of STD RSF table.
- Further steps:
- Step 4 pulsed mode with separate RSF table. First steps done for 33 elements.
- Step 5 Round Robin study to create new, averaged STD RSF table in the future. Materials e.g. Ni alloys.
- Step 6 rename STD RSF into e.g. "Generic RSF".



Vielen Dank für die Aufmerksamkeit!



joachim.hinrichs@thermofisher.com

