

Optical emission spectra of metallic nitrides in the Grimm- type glow discharge

Petr Šmíd^a, Zdeněk Weiss^b and Edward B. M. Steers^c

^a *Imperial College London, UK*

^b *LECO Instrumente Plzeň, Czech Republic*

^c *London Metropolitan University, UK*

Copyright:

Dieser Vortrag darf nur nach Absprache den Autoren zitiert werden

Goals of this paper:


- to compare optical emission spectra of metallic nitrides with corresponding metals and to test validity of matrix independent emission yield concept
- to compare optical emission spectra of the discharge when:
 - nitrogen enters the discharge from a metallic nitride
 - nitrogen is introduced into the discharge as a trace molecular gas

Experimental – metallic nitrides

- commercial instrument LECO GDS500A with a Grimm-type glow discharge sources, 4 mm diameter anode tube
- spectrometer equipped with CCD detector arrays, wavelength range from 165 to 460 nm and spectral resolution of 70 pm (310–460 nm), resp. 85 pm (165–310 nm)
- spectra of metallic nitrides – nitrided steel, CrN, TiN and ZrN (stoichiometry checked with RBS) – compared with spectra of corresponding metals – steel, Cr, Ti and Zr
- constant dc voltage (700 V) and constant current (20 mA)

Experimental – metallic nitrides

Qualitative depth profile of the nitrided steel sample




Experimental – Ar/N₂ discharges


- Grimm-type glow discharge sources, 8 mm diameter anode tube
- Fourier Transform Spectrometer with PMT detection, wide wavelength range (135–900 nm) and high resolution (e.g. 0.85 pm at 300 nm)
- Fe and Ti samples sputtered in the glow discharge running in pure Ar and Ar with small amounts of N₂ (0–2 % vol./vol.)
- constant dc voltage (700 V) and current (40 mA)

Results – metallic nitrides





Spectrum of a steel sample (DIN X38-CrMoV51)









Spectrum of a **nitrided** steel sample (DIN X38-CrMoV51)



Results – metallic nitrides



Results – metallic nitrides vs. Ar/N₂ discharges



Results – argon lines

Intensity of an emission line of a wavelength λ :

$$I_{\lambda} \propto n_2 A_{21} \frac{h\lambda}{c}$$



Plasma gas emission line:

$$I_{Ar,\lambda} \propto p \cdot f_{exc}(\lambda, p) \cdot A_{21} \frac{hc}{\lambda}$$



Intensity ratio:

$$\frac{I_{Ar,\lambda}^{MN}}{I_{Ar,\lambda}^M} \propto \frac{p_{MN}}{p_M} \cdot \frac{f_{exc}^{MN}(\lambda, p_{MN})}{f_{exc}^M(\lambda, p_M)}$$

Results – argon lines

$$\frac{I_{Ar,\lambda}^{MN}}{I_{Ar,\lambda}^M} \propto \frac{p_{MN}}{p_M} \cdot \frac{f_{exc}^{MN}(\lambda, p_{MN})}{f_{exc}^M(\lambda, p_M)}$$

Comparison of average argon intensity ratios with gas pressure ratios:

Metallic nitrides:

Samples	Intensity ratios of argon emission lines			Gas pressure
	Mean	STD	Number	ratio p_{MN}/p_M
Nitrided steel/steel	0.96	0.03	60	1.04
CrN/Cr	0.80	0.05	60	1.05
TiN/Ti	0.78	0.05	60	1.19
ZrN/Zr	0.80	0.07	60	1.04

Ar+N₂ mixtures:

	% N ₂ vol./vol.	Intensity ratios of argon emission lines			Pressure
		Mean	STD	Number	ratio
Fe	0.05	0.97	0.05	60	1.00
	0.18	0.92	0.18	60	1.02
Ti	0.07	0.98	0.06	60	1.06
	0.11	0.92	0.12	60	1.15

Results – sample metal lines

Intensity of an emission line of a wavelength λ :

$$I_{\lambda} \propto n_2 A_{21} \frac{h\lambda}{c}$$



Sample emission line:

$$I_{E,\lambda} = R_{E,\lambda} c_{E,M} q_M$$



Intensity ratio:

$$\frac{I_{E,\lambda}^{MN}}{I_{E,\lambda}^M} = \frac{R_{E,\lambda}^{MN}}{R_{E,\lambda}^M} \cdot \frac{c_{E,MN} q_{MN}}{c_{E,M} q_M}$$

Results – sample metal lines (metallic nitrides)

$$\frac{I_{E,\lambda}^{MN}}{I_{E,\lambda}^M} = \frac{R_{E,\lambda}^{MN}}{R_{E,\lambda}^M} \cdot \frac{c_{E,MN} q_{MN}}{c_{E,M} q_M}$$

For Matrix Independent Emission Yield:

$$R_{E,\lambda}^{MN} = R_{E,\lambda}^M \Rightarrow \frac{I_{E,\lambda}^{MN}}{I_{E,\lambda}^M} = \frac{c_{E,MN} q_{MN}}{c_{E,M} q_M}$$



Average intensity ratio vs. flux ratio

	Intensity ratio	Min	Max	STD	No. lines	Flux ratio	STD
Nitrided steel/steel	0.84	0.71	0.95	0.05	134	---	---
CrN/Cr	0.42	0.35	0.55	0.05	46	0.39	0.05
TiN/Ti	0.47	0.34	0.66	0.05	173	0.49	0.11
ZrN/Zr	0.41	0.31	0.61	0.05	239	0.43	0.09





Results – sample metal lines (metallic nitrides)

For Matrix Independent Emission Yield:

$$R_{E,\lambda}^{MN} = R_{E,\lambda}^M \Rightarrow \frac{I_{E,\lambda}^{MN}}{I_{E,\lambda}^M} = \frac{c_{E,MN} q_{MN}}{c_{E,M} q_M} \neq f(\lambda)$$



Results – sample lines (metallic nitrides vs. Ar/N₂ discharges)




Results — nitrogen lines (metallic nitrides vs. Ar/N₂ discharges)


Emission yields for nitrogen atomic lines:

	$R_{\text{NI } 174.273}$	STD	$R_{\text{NI } 174.525}$	STD
<i>Nitrided steel</i>	---	---	---	---
<i>CrN</i>	3.9	0.4	1.8	0.2
<i>TiN</i>	4.1	0.9	1.7	0.4
<i>ZrN</i>	3.5	0.5	1.6	0.2

Intensity of N I 174.525 nm vs. N₂ concentration in Ar+N₂:




Results – nitrogen molecular bands (metallic nitrides)



Results – nitrogen molecular bands (Ar/N₂ discharges)

Intensity of N₂ 337.13 nm molecular band vs. N₂ concentration in Ar:



Acknowledgement

Petr Šmíd would like to thank to

- the Royal Society for a one year RS/NATO postdoctoral fellowship at Imperial College London
- the Association of British Spectroscopists Trust for funding his attendance at this meeting